## **Leuschner Operational Instructions – Quick Start - revised 2-23-2023**

# **Beginning of Night Procedure**

#### LOG in to WEBCAM

- 1. Chrome web browser: search for Amcrest Webview app [EOL -unsupported, but required]
- 2. Amcrest Webview: search leuschner.berkeley.edu:82
- 3. Log in window: User: talon Pswd: {see Admin} [Under Settings: Turn on IR light if OFF]

#### LOG In to TALON

- 1. Terminal: ssh -X (or -Y) talon@leuschner.berkeley.edu
- 2. Terminal: ssh -X (or -Y) talon@192.168.1.195
- 3. Terminal: firefox [opens firefox web browser in XQuartz]
- 4. Firefox: Go to 192.168.1.210:81 [opened talon log in window within firefox]
- 5. DDL Controller: Click on SBIG CCD power to ON, then logout
- 6. xQuartz: Close firefox window in XQuartz.
- 7. Terminal: *killTel -all*, then *StartTel -all* [opens all xobservatory windows plus flip mirror window (camera program opened)]

### HOME ALL TELESCOPE FUNCTIONS

- 1 xObservatory: Click on Find Home
- 2 xObservatory: **Home** window: Click on **Scope Dec/Alt** [Telescope finds home position-Confirm in Webcam and Text window]
- 3 **xObservatory**: Click on **Paddle**: Press Y button until **Dec** = **+15°** [Confirm in webcam and in xObs Current position]
- 4 xObservatory: Home window: Click on **Scope RA/Az** [Telescope finds home position-Confirm in Webcam and Text window]
- 5 **xObservatory**: Click on **STOW**. **xObs Query window YES** [Telescope finds STOW position-Confirm in Webcam and Text window]
- 6 **xObservatory:** Home window: Click on **Dome** {Dome rotates to Home position Confirm in webcam and text window]
- 7 xObservatory: Home window: Click on Focus [Secondary mirror moves to home position Confirm in webcam, xObs, and text field 25,600]

#### INITIATE CAMERA FILTER PROCESS

- 1. Terminal: ps aux/grep camera [ outputs process numbers for camera ]
- 2 Terminal: *kill {camera process#}* [camera xwindow closes]
- 3 Terminal: cd /home/talon/xmccd-3.1.0/xmccd1 [change directory]
- 4 Terminal: -./xmccd1 [error message followed by opening camera filter window]
- 5 Xmccd1: set camera **filter**, {select 1=1, 2=R, V=V(green), 4=B, or 5=U},
- 6 XQuartz: close xmccd1 window
- 7 Terminal: cd /home/talon
- 8 Terminal: camera & [camera program re-opens]
- 9 Camera 3.44 Xwindow- set Exposure : Dropdown: Expose: Choose Setup...
- 10 Exposure Setup Xwindow: Binning = 2:2; Duration = {exposure time} Second(s) [Press Start to take a CCD image after readout, image appears in Zoom window]

#### **BEGIN OBSERVATIONS**

### Select Instrument in Flip Mirror:

Port 1: SBIG CCD Camera; [See Leuschner Operational Manual for details]

Port 2: Visual Eyepiece Focuser [Onsite Use only];

Port 3: SFSU eShel Spectrograph – Requires separate LOG IN Procedure simultaneously.

[See <u>SFSU Spectrograph Operations Manual</u> for details]

## Open Dome and Mirror Cover:

- 1 xObservatory: Roof Open [Dome shutters open and Mirror Cover open Confirm in webcam and Text window] -NOTE: Dome shutters open/close only if Dome is in Home position for power, but Mirror covers will always open/close regardless of Dome position.
- 2 xObservatory: Roof: Check Auto box xObs Query window YES: Status indicator turns green from red [Move the open shutters to the telescope position automatically and moves the dome as telescope tracks a target – Confirm in webcam]

## **End of Night Procedure**

#### STOW all TELESCOPE FUNCTIONS

- 1 **xObservatory:** Click on **STOW. xObs Query window YES** [Telescope finds STOW position-Confirm in Webcam and Text window]
- 2 **xObservatory:** Home window: Click on **Dome** {Dome rotates to Home position Confirm in webcam and text window}
- 3 xObservatory: Click on Roof: Close {Dome Shutters close; Mirror covers close Confirm in webcam and text window]

#### SHUTDOWN all PROGRAMS and LOGOUT of TALON

- 1 Terminal: *firefox* [opened firefox within XQuartz]
- 2 XQuartz: Firefox: Go to 192.168.1.210:81 [open talon log in window]
- 3 Firefox: User: talon PSWD {See Admin}
- 4 Firefox: DLL Controller: Click on SBIG CCD power to **OFF**, then **logout**
- 5 XQuartz: Close *firefox* window.
- 6 Terminal: killTel -all
- 7 Terminal: **exit** <u>talon@tcs-desktop</u> -> **exit** <u>talon@leuschgate</u>
- 8 Chrome: logout of Amcrest Webview Close Chrome